

## GM-2Z GEM ZOOM STEREO MICROSCOPE

### composed of:

- EMZ-2** Zoom stereo microscope body with zoom objectives 0.7X to 4.5X. Magnification from 7X to 90X (using SWF10X and SWF20X eyepieces)
- ABZH/1** Stand with transmitted illuminator, 6V 20W Halogen lamp and overhead illuminator, 4W fluorescent lamp
- MA502** Super widefield 10X eyepiece, FN 23mm (paired)
- MA504** Super widefield 10X eyepiece, FN 11.5mm (paired)
- MA555** Gem clamp
- MA556** Darkfield attachment
- MA557/01** Iris diaphragm with adapter

## GM-2T GEM STEREO MICROSCOPE

### composed of:

- EMT-2** Stereo microscope body with 1X and 3X objectives in turret changer. Magnification: 10X, 20X, 30X and 60X (using SWF10X and SWF20X eyepieces)
- ABH/1** Stand with transmitted illuminator, 6V 20W Halogen lamp and overhead illuminator, 4W fluorescent lamp
- MA502** Super widefield 10X eyepiece, FN 23mm (paired)
- MA504** Super widefield 20X eyepiece, FN 11.5mm (paired)
- MA555** Gem clamp
- MA556** Darkfield attachment
- MA557/01** Iris diaphragm with adapter

## GM-3T GEM STEREO MICROSCOPE

### composed of:

- EMT-3** Stereo microscope body with 2X and 4X objectives in turret changer. Magnification: 20X, 40X and 80X (using SWF10X and SWF20X eyepieces)
- ABH/1** Stand with transmitted illuminator, 6V 20W Halogen lamp and overhead illuminator, 4W fluorescent lamp
- MA502** Super widefield 10X eyepiece, FN 23mm (paired)
- MA504** Super widefield 20X eyepiece, FN 11.5mm (paired)
- MA555** Gem clamp
- MA556** Darkfield attachment
- MA557/01** Iris diaphragm with adapter

## GM-4T GEM STEREO MICROSCOPE

### composed of:

- EMT-4** Stereo microscope body with 1X and 4X objectives in turret changer. Magnification: 10X, 20X, 40X and 80X (using SWF10X and SWF20X eyepieces)
- ABH/1** Stand with transmitted illuminator, 6V 20W Halogen lamp and overhead illuminator, 4W fluorescent lamp
- MA502** Super widefield 10X eyepiece, FN 23mm (paired)
- MA504** Super widefield 20X eyepiece, FN 11.5mm (paired)
- MA555** Gem clamp
- MA556** Darkfield attachment
- MA557/01** Iris diaphragm with adapter

### SPARE BULBS:

- MA309/05** Fluorescent tube, 4W
- MA311** Halogen lamp, 6V 20W



GM-2Z



GM-T Series

**ISSCO**

Ph.(02) 9743 2588  
Fax.(02) 9736 2835

SIZE & WEIGHT : 170mm(W)x290mm(D)x410mm(H),6.3kg



# GS2000

SERIES  
GEM STEREO  
MICROSCOPES

## ECONOMICAL GEM STEREO MICROSCOPE WITH DARKFIELD / BRIGHTFIELD BAFFLE

### STANDARD FEATURES

- \* High resolution optics that deliver sharp, clear images and bright accurate colors
- \* Turret-type, two position magnification changer
- \* Darkfield/brightfield baffle
- \* Binocular body inclined at 45° for ease of observation
- \* Widefield 10X eyepieces with eyeguards
- \* Interpupillary adjustment between 51mm and 75mm
- \* Dioptic adjustment of ±5mm on one ocular tube
- \* Objectives and magnification according to models
- \* Pole type stand with heavy metal base, with focusing block
- \* Built-in incident light illuminator 12V 10W
- \* Transmitted daylight illuminator direct from mains, 115V or 220/240V 20W
- \* Separate ON/OFF switches for incident and transmitted illuminators
- \* Gem clamp
- \* Supplied in styrofoam case with dustcover and manual

### ORDERING INFORMATION

MODEL	OBJECTIVES	EYEPIECES	MAGNIFICATION
GS2100	1X, 2X	WF10X FN20	10X, 20X
GS2200	1X, 3X	WF10X FN20	10X, 30X
GS2300	2X, 4X	WF10X FN20	20X, 40X

SIZE & WEIGHT: 115mm(W)x170mm(D)x270mm(H), 2.8kg

### OPTIONAL ACCESSORIES

- MA554 Iris diaphragm with adapter
- MA102 Widefield eyepiece 10X, FN 20 (standard package)
- MA103 Widefield eyepiece 15X, FN 15
- MA104 Widefield eyepiece 20X, FN 13
- MA555/05 Gem clamp (included)
- MA589/100 Spare lamp, 115V 20W Daylight
- MA589/200 Spare lamp, 220/240V 20W Daylight
- MA111 Spare lamp, 12V 10W, Tubular bulb (incident light)
- MA505 Eyeguards (paired)



GS-2000 Series

Agent :

**ISSCO**

Ph.(02) 9743 2588  
Fax.(02) 9736 2835



**MEIJI TECHNO CO.,LTD.**

6, OI-670, OI-machi, Inuma-gun  
Saitama 356, Japan  
Phone : 492-67-0911  
Fax : 492-69-0691, 492-69-0692  
Telex : 25573 MEIJTCHN  
Cable : MEIJITECHNO OI

**MEIJI TECHNO AMERICA**

2186 Bering Drive  
San Jose, CA., 95131, USA  
Phone : 408-428-9654 Fax : 408-428-0472  
Toll free : 800-832-0060